



Order Form - Comprehensive AFO

Date: _____ P.O. # _____
 Account Name: _____
 Billing Address: _____

 City: _____
 State: _____ Zip: _____
 Phone: _____ Fax: _____
 Email Address: _____

Ship To (if different than billing address):
 Company Name: _____
 Address: _____
 City: _____ State: _____ Zip: _____
 Phone: _____ Fax: _____

Ship Method (charges may vary):
 Next Day 2 Day 3 DAY Ground (std.)

Standard Specifications:

- 3/16" Black or Natural Polypropylene unless otherwise specified
- Color:** Black Natural
- Brace:** Lined Unlined
- 1 Strap at top Add instep strap
- Standard - non-articulated (as shown above)
- Rigid Semi-Rigid Flexible (pt. wt. required)*

Plastic Length: To mets (Std.) To sulcus End of toes
Top Cover length: To sulcus End of toes

*Options:

- *Tamarack Dorsi Assist (Pre-Flexed 75 durometer - Std.)
- *Dorsi assist with plantar flexion stop
- *Limited range of motion control (straight Tamarack hinge)
- *Additional charges may apply

Cast Mold Information:

A cast mold of the patient's leg should be taken with a "Bermuda Casting Sock" in a non-weight bearing position with the foot and calf relationship at 90°. The patient's foot can be placed on a soft foam cushion after its cast. Position the foot and ankle in the desired position of the finished orthosis. AFO cast mold should go up to 1-2 inches below the head of fibula. A foot board may be used if desired (NOTE: Please NO Webriil/cotton casting material!)

Important Notes:

- Indicate all bony landmarks with indelible pencil
- All measurements must be taken for reference
- An above calf cast MUST be taken for these orthoses.

With each shipment, we provide an STS casting sock free of charge.

- Bermuda or Mid-calf
- Medium Large X-Large None at this time

Patient Info:

Patient Name: _____
 Age: _____ Sex: M F Ht: _____ Wt: * _____
 Occupation: _____
 Primary Activities: _____

 Diagnosis: _____
 Shoe Size: _____ Type of Shoes: _____

Brace Side:

Left Right Bilateral

Special Instructions:

Note: Above calf cast is ABSOLUTELY necessary.

Mark Accommodation Location(s) on Cast.

